

Lab Notebook

Adam Schonewille 31590136

Akshiv Bansal 35072131


Day 1 2018 03 06

Optical Set-up, Magnification Observations and Wire Mesh Experiments:

Optical Setup and Alignment

Mesh that we used has 40 lines/cm

Fourier Optics Lab



Moving the collimator lens does not effect magnification or focus. It does seem to vary the intensity on the outer edges of the circle.

1. Why is the spatial filter important for this lab?

The primary purpose of the spatial filter is to remove any higher order spatial modes that maybe present in the laser output. This helps to create a single Gaussian beam, which behaves more predictably and is easier to model.

2. What is the best object to use in order to define an object and image plane combination? Why? Use this object to find an object plane and the corresponding image plane.

First we used the pointed rod to align the spatial filter and laser combination, ensuring our beam was parallel to the rail at all points. Then to find the object plane we placed the razor blade with the clear plastic wedge attached to in approximate location of the object plane. Then we made slight adjustments until we could produce a sharp image on the screen.

3. What is the best object to use in order to find the Fourier transform plane? Use this object to find the Fourier transform plane

To find the Fourier transform plane, we used the mesh filter because it produces a very identifiable Fourier transform. With this object in the object plane we moved a screen into the approximate location of the Fourier transform plane and adjusted it until it was sharp. We then recorded the location of both of these planes.

We placed the spatial filter at the following settings:

Image plane

147.5+358.2+4.2

1.3-167

Absolute means from the start of measurement rail, offsets are already included in reported distances. Raw measurements are taken from the back of the stage and the offsets correct to the middle of the optical apparatus.

For ease of reconstruction:

Object	Set-up Distance at back of rail slider (cm)
Spatial Filter	33.3
Transform Lens	167
Mesh	108.5
Collimating Lens	55.8

Object	Set-up Distance at back of rail slider (cm)
Imaging Lens	265.5
Object Plane from knife edge	107.1

The distances from the back of several objects to the center of the lenses/apertures is recorded for relative measurements. This was the set-up distances can be related to the locations of action and used in the calculations later on.

Object	Offset Distances (cm)
Spatial Filter Offset (Δ)	11.5
Transform Lens back offset (Δ)	1.3
Mesh Offset (Δ)	3.2
Collimating Lens Offset (Δ)	3.05
Imaging Lens Offset (Δ)	3.5

Doing all these calculations gives the true locations of all the components in this optical setup:


Object	Absolute Distance for calculations (cm)
Spatial Filter	44.8 ± 0.5
Transform Lens	168.3 ± 0.5
Image Plane	509.9 ± 0.5
Mesh	111.7 ± 0.1
Collimating Lens	58.85 ± 0.1
Object Plane from knife edge	107.1

Imaging lens was placed after the Fourier transform plane

Mesh Filtering Experiment

Grid spacing from paper: 1/8"

Mesh that we used has 40 lines/cm



1. From the number of wires per cm given on the mesh aperture, and the spacing of the image wires on the screen, estimate the magnification of this system.

This data was retaken later in the lab and is mentioned later in this book.

3. Separately, use an additional lens to generate a magnified image of the Fourier transform plane onto the screen and take a picture.



- . What do the bright spots in the Fourier transform plane represent?

Each point in the pattern represents a unique spatial frequency, and is therefore the focus of all parallel rays in the object space making a certain angle with the optic axis

Other questions located later in the lab where this part was redone


9. Rotate the mesh around in the mount and lift it up and down. How does moving the mesh change the magnified image and the Fourier transform? Why do some types of movement change the image but appear not to change the Fourier transform of the image? What might be happening?

For the object in the image plane the image translates and rotates in the same manner as the actual movement of the object in the object plane. However, the Fourier transform of the wire mesh dims in areas for some fluctuations of the object in the object plane. The absence of some changes may be due to symmetries in the wire mesh when it is Fourier transformed.

4. What do the bright spots in the Fourier transform plane represent?

The

5. Rotate the mesh around in the mount and lift it up and down. How does moving the mesh change the magnified image and the Fourier transform? Why do some types of movement change the image but appear not to change the Fourier transform of the image? What might be happening?




10


6. Using the mesh object, how can you produce only vertical lines? Explain how and why this is happening. Provide a picture.

ToFourier transform image, calculate the spatial frequencies that are present. (Please be careful with units.)


7. What wire spacing do these spatial frequencies correspond to this correspond to?
8. How does this we can use a vertical slit aperture compare to what you expect?
9. Using the Fourier transform plane to only allow vertical components of the image through. This will produce a grid with only vertical components
mesh object, how can you produce only vertical lines? Explain how and why this is happening. Provide a picture.



Save Figure




Looking at the Fourier transform here makes it clear that only these vertical components survived. This is made more obvious if we look at the 3D transform.



10. Try producing lines at 45° and horizontal. If you are able to produce these lines, explain how you did it and explanation of why its works. Include a picture with your explanation.







Dark-field Image

A Dark-field Image can be created by using a Fourier transforming lens and placing a point aperture right at the Fourier transform plane (slightly blocking the focus) and viewing the resulting image in the image plane.

Normal Image



Dark Field Image




1. What is meant by a “dark-field image”? What aperture can you use to produce a dark field image?

The dark field image is missing the central spatial frequency, in the Fourier transform sense. This means that the zeroth order spot in the Fourier plane is removed. We can produce a dark field image by using the aperture with a dark spot in its center and passing the beam focused to a point through that spot.

2. Which object will give the best visibility of the “dark-field image” effect? Justify your answer

In our opinion the razor blade with a plastic wedge glued to it provides the best dark field image. This is because the dark field illuminates triangle shape, which is invisible in the normal image. The dark field image clearly shows the outline of the triangle, which was absent in the first picture.

We also tried dark field imaging with a mesh transform. As expected it removed the central bright spot in the diffraction pattern.



Day 2 2018 03 13


Character Recognition, Phase-Contrasting and Diffraction:

Character Recognition:

In today's experiments we placed a sheet of grid paper onto the image plane to be able to see the difference in skew that the camera has on the image plane. This is because the camera views the image plane at an angle and so the bottom of the field of view is smaller than the top. First we set up the optical rail as follows to view the NOZON object in the image plane:

- **Experimental set-up:**
- Collimating Lens
- NOZON aperture in object plane
- Fourier Transform Lens

We were able to get a clear image of the NOZON aperture:




We can also invert the image to get a clearer image:

NOZON
NNOZN
NONON
NZONN
NOZON


A white rectangular image containing clear black text. The text is identical to the one in the previous image but is much more legible due to the inversion. It consists of five lines of text: 'NOZON', 'NNOZN', 'NONON', 'NZONN', and 'NOZON'.

We attempted to adjust the gain so that the image can be seen against the grid of paper to view the camera skew:




Next we added the Imaging lens into the set up to view the fourier transform:

- **Experimental set-up:**
- Collimating Lens
- NOZON aperture in object plane
- Fourier Transform Lens
- Imaging Lens to see the FT of the NOZON object




Can also invert the image to save on ink and see more clearly what is going on:




The FT of the NOZON aperture looks similar to an asterisk (*) and in order to



Next we put the Lazy X in




Again we invert the image:



Because the NOZON aperture is much less regular than the mesh filter the Fourier transform image is significantly more complex. We need many more spatial frequencies to achieve the sharp corners involved in the NOZON filter. Looking at the Fourier transform image we see that these higher order frequencies and peaks represent the sharp edges of the image.



The reason we can use the lazy x to achieve the correct spatial filtering is we want block the smoother frequencies to get rid of the "O"s and the flat bars in the "Z"s. We we try to maintain the peaks and this is possible by applying a filter which removes just those frequencies.




Looking at the lazy x filtered image, the “round” frequencies have been attenuated into the noise, while the sharp frequencies maintain their relative amplitude.

Phase-Contrast


Experimental set up:

Object	Set-up Distance at back of rail slider (cm)
Spatial Filter	33.3
Collimating Lens	102.3
Phase Grating	148.4
Transform Lens	204.0
Wedge Knife Blade	236.2
Imaging Lens	248.5

Next we collected Images for the phase grating with and without the FT lens and with and without the "Phase Plate" (vertical razor blade with a wedge of clear plastic taped to it into the Fourier Transform plane). These images are shown below.




Phase grating with Fourier Transform, but not Phase Plate. Here the image looks homogeneous and the difference in phases are not noticeable.



By moving the Phase Plate into the Fourier transform field we see that the differences in phase become clearer as they are distinguished. Here we expect to see a grid, but something is wrong with our phase grating. We talk more about this below.

Without the Fourier transform lens but with the Phase Plate we should be able to see the vertical lines of the difference in phase on the image plane caused by the phase grating. The following image shows these results.



This image shows the phase difference is distinguishable and we were able to measure the difference in phase.

Bottom-right fringes separation (between successive bright nodes):

$$5 \pm 0.5 \text{ mm}$$

It should be noted that the phase grating that we were using had visible problems as these images show the image is skewed. There should be completely vertical lines in the phase grating image on the screen but instead they were curved. We asked the TA's about this and the consensus was that the phase grating must be either damaged or dirty. Because of this, we could not get great data and decided to not make this the primary focus of this lab.

Bottom-right fringes separation (between successive bright nodes):

$$5\text{mm}$$

Top-middle-right fringes separation (between successive bright nodes):

$$3.5\text{-}4\text{mm}$$

1. Explain how the grating causes a spatially dependent phase in the original light wave.

Phase contrast imaging is a technique by which spatial filtering is used to transform the phase variation across the wave into amplitude variations which can then be detected. We can say that the phase (shift) of light that passes through a medium is given by:

$$f(x, y) = e^{i\theta}$$

Which is approximated by:

$$f(x, y) \approx 1 + i\theta$$

Then the intensity is:

$$|f(x, y)|^2 \approx 1 + i\theta - i\theta - i^2\theta^2 = 1 + \theta^2$$

Visually, us, cameras and many other devices detect light by direct measurement of the light intensity. And since $\theta \ll 1$ the phase difference is not usually distinguishable.

If we were to phase shift the DC component by 90° (by a factor of i) we should be able to distinguish the differences in phase much easier. Here $f(x, y)$ becomes:

$$f(x, y) \approx i + i\theta$$

$$|f(x, y)|^2 \approx -i^2 - i\theta - i\theta - i^2\theta^2 = 1 + 2\theta + \theta^2$$

Even though 2θ is still much smaller than 1, it gives enough contrast to be able to distinguish the phase difference.

2. Insert the vertical razor blade with a wedge of clear plastic taped to it into the Fourier Transform plane. How does this “phase plate” allow you to modify the phase of the DC component while leaving all (most) of the other components unaffected?

Since the DC term has no spatial variation, representing just a flat and featureless illumination from the plane wave, it will be focused to a point in the Fourier Transform plane and comprise the zeroth order of the diffraction pattern. By simply inserting a thin dielectric layer (film) at the central spot in the transform plane the phase of the uniform illumination will be altered. A quarter wavelength phase change would produce the best results (90° DC phase shift).


3. Using the phase-contrast image of the phase grating, compute the magnification of your imaging system. How does this compare to your previous measurement with the mesh?

If the grating has about 0.3 lines/mm and we observed that there was about 5 mm/line that gives a magnification of about 1.5


This value for magnification is 4.32 times smaller than the magnification found for the mesh experiment. This is possible however because we were not using the same optical setup for these two days so the positioning of the optical components will have varied by some amount.

This method can be used to determine if there is varying thickness in a media, but not able to determine the exact thickness of the refractive medium. This is because a phase shift of 2π will not show up after the phase filtering is applied. Relative error will however be able to be used to distinguish changes in thickness as the length and difference in refractive index will phase shift the light by some amount and will become visible after phase filtering in the FT plane. Experiment to find Magnification


To measure the magnification of the image, we used graph paper as a reference and plotted the intensity of the light obtained from the CCD camera image.



Here is the slices that were taken for processing data



This is the data for both slices. In blue is the reference slice for pixels to grid spacing measurements. The dips in this signal are the black lines seen on the right of the image. In red is the brightness plot of the wire mesh image. the peaks are the dots seen from the screen.



Redo of Mesh Filtering Experiment to find Magnification

The grid paper that we were using had spacing that could be measured via image processing and compared to the measured value (with a ruler) to obtain a relationship between the camera and the actual spacing.

Grid spacing from paper: $1/8"$ \rightarrow 3.175 mm

Mesh that we used has 40 lines/cm \rightarrow 4 lines/mm

This mesh can also be described as $250 \mu\text{m}/\text{line}$ (spacing between lines)

Peak Pixel Location	Difference Between Peaks	Mean (Pixels per Peak)
408	N/A	
414	6	
421	7	
426	5	
433	7	
440	7	
446	6	6.33

This data shows that there is a wire/line every 6.33 pixels on the screen

Dip Pixel Location	Difference Between Peaks	Mean (Pixels per Dip)
436	N/A	
448	12	
460	12	
472	12	
485	13	
498	13	12.4

The reference image gives the spacing relation between mm and pixels

There is 3.175 mm/dip and there is also 12.4 pixels/dip

Therefore there is about 3.9 pixels/mm

From this we can assume that on the screen there is 1.62 mm/line


This is a magnification by a factor of 6.487 times larger.

We didn't measure the exact lengths between certain optical components so a lens analysis of this number is difficult however we do have some numbers to see if our result is in the right ballpark. If the Image plane is at 509.9 cm, the Fourier transform lens is at 168.3 cm and the object plane at 111.7 cm we can determine that the magnification is roughly:


$$M = (509.9 - 168.3) / (168.3 - 111.7) = 6.04 \pm 0.29$$

Where the error in the location of each plane is roughly $\pm 5\text{ cm}$. Here we see that the magnification is very close to the magnification that we got before.


7. Using the Fourier transform image, calculate the spatial frequencies that are present. (Please be careful with units.)




This is the Image of the FT that was analyzed. There are reference bars to see the spacing in this field.



These pixel intensities were plotted and shown above. The blue curve shows the intensity of the pixels going vertically through the reference lines (dips) and then through 3 spots of the FT (peaks). The red curve is the pixel intensity vertically through the center of the image. These lines have a lot of noise and these peaks/data was ignored. The green curve shows the horizontal pixel intensities through the center of the image horizontally. The left-most peaks were observed alone.



Horizontal spacing of the FT image of the wire mesh.



Reference spacing measurements for this image.

Peak Pixel Location	Difference Between Peaks	Mean (Pixels per Peak)
68	N/A	
110	42	
153	43	
195	42	
238	43	42.5
Peak Pixel Location	Difference Between Peaks	Mean (Pixels per Peak)
152	N/A	
192	40	
234	42	41

This data shows that there is a peak every 42.5 pixels on the screen

Dip Pixel Location	Difference Between Peaks	Mean (Pixels per Dip)
42	N/A	
52	10	
63	11	
75	12	
86	11	
97	11	11

The reference image gives the spacing relation between mm and pixels

There is 3.175 mm/dip and there is also 11 pixels/dip

Therefore there is about 3.46 pixels/mm

From this we can assume that on the screen there is 12.28 mm/peak

using the magnification we found before, we can determine that at the wire mesh this corresponds to a distance of 1.89 mm/peak at the object plane.

Day 3 2018 03 20

Fraunhofer and Fresnel Diffraction Regimes:

Fresnel Magnification

We had the idea to take a dark-field image of the setup with the blades. Previous Dark-field Images have produced only the

Aperture settings:

Blade aperture separation: 21.20 ± 0.01 mm (during measurement)

Blade aperture separation: 23.40 ± 0.01 mm (fully closed)

Set Aperture

Measured distance with a ruler:

On paper 12.7 ± 0.2 mm

On screen 13.0 ± 0.2 mm

Measured with Image Processing:

For all Optical Components

For no Optical Components:

Measured slit width on image plane: 5.5 ± 0.5 mm

Set Blade Aperture Position: 20.15 ± 0.5 mm

Blade aperture separation: 23.40 ± 0.01 mm (fully closed)

The screen in this setup has a slit-screen distance of 3.85 ± 0.01 m

We remove all lens from the optical set up and leave only a light source, a variable slit aperture, an imaging screen and the CCD.

We develop a theoretical diffraction pattern, to predict a diffraction pattern seen at a screen in the far field(Fraunhofer) and near field(Fresnel) limits. While the mathematical simplicity of the Fraunhofer case is tempting, we simply develop the Fresnel patterns, and allow our limiting values to determine the pattern in the Fraunhofer case without any additional simplification.

We start by defining a quantity Δv , which is a dimensional parameter representing the “nearness” of a diffraction pattern. If Δv is large we are squarely in the Fresnel limit, while if we are in the situation where $\Delta v \ll 1$ we are in the Fraunhofer.

$$\Delta v = w \sqrt{\frac{2}{R\lambda}}, \text{ where } w = \text{slit width}, R = \text{slit to screen distance}, \lambda = \text{operating wavelength}$$

Then using the Fresnel Reflection equations and letting C be some constant of proportionality we get:

$$I(z) = C \int_{v_1}^{v_2} \cos \frac{\pi x^2}{2} dx + C \int_{v_1}^{v_2} \sin \frac{\pi x^2}{2} dx, \text{ where } v_1 = -(z + 0.5)\Delta v, v_2 = -(z - 0.5)\Delta v$$


As expected in the far-field where $\Delta v \ll 1$ we recover:



$$I(z) = C(\Delta v)^2 \left[\frac{\sin \left(\frac{z\pi(\Delta v)^2}{2} \right)}{\frac{z\pi(\Delta v)^2}{2}} \right]^2 = C(\Delta v)^2 \frac{\sin^2 \beta}{\beta^2}$$



To see the effect of change Δv we can either move the screen back and forth while holding everything else constant, or we can adjust the width of the slit (I suppose we could also change the operating wavelength, but that is not possible in a practical sense). Since we want to maintain the imaging setup we have for screen, we choose to modify the slit width in order to adjust the diffraction regime. We start by defining the slit's closed point where no light from the laser is visible, and we count this as the background intensity. We will use this reference for the differential slit width, and this image intensity as the background subtraction to compensate for background optical noise. Then we slowly open the slit and take images on the CCD, at various points. One could take a regular series of slit width data, but since



interesting data is only really found when we squarely in the Fresnel or Fraunhofer regions, we elected to images and slit width when we saw interesting fringes. At minimum we need one image in the far field, one in the near field, and lastly one at the transition point.


Using the MATLAB script, which you can find here: <https://github.com/akshivbansal/phys408FourierOptics/tree/master/processedData/Diffraction> we averaged 10 pixels in the center of each image. We then compared these diffraction patterns with their Fourier transforms and their predicted diffraction patterns.





















Diffraction Pattern in Transition $\Delta v = 3.16$







Diffraction Pattern in Transition $\Delta v = 2.52$







Diffraction Pattern in Transition $\Delta v = 1.78$







Diffraction Pattern in Transition $\Delta v = 1.1$







Diffraction Pattern in Transition $\Delta v = 1.1$







Diffraction Pattern in Transition $\Delta v = 0.68$






Diffraction Pattern in Transition $\Delta v = 0.43$





Diffraction Pattern in Transition $\Delta v = 0.32$





In general we see that that the data is better modeled as we approach the far-field limit. This is because the inherently discrete nature of our model and the approximations made in the model matter less as we move to the far field limit and recover a cardinal sine function. Importantly, the theory does a good job of predicting the shape and relative magnitude of the intensity of the diffraction pattern. While we can't use the theoretical values to perfectly model the diffraction in these cases we are still able to get a good sense for the intensity profile in most situations. All of our plots suffer from horizontal scaling issues. This cause by a difference in the imaging screen and the actual size of the image produced on it. There is also a magnification correction we applied, which has error associated with it, this combined with the skew of the image , likely causes difference between real and imaged pixel size. Over all, the theory of Fresnel diffraction seems to good job of predicting the observations made in the lab.